[B] INORGANIC QUALITATIVE ANALYSIS

Inorganic Qualitative analysis deals with the detection of an unknown powder or solution by the systematic chemical methods. The substance under analysis is generally composed of two parts. Each one is known as radical. Like, positive radical or cation (+Ve) and negative radical or anion (-Ve).

Usually the substance under analysis is subjected to:

- (A) Dry tests (For substance in the solid state only).
- (B) Wet tests (For substance in the aqueous state only and for the solution prepared from solid salts).

General procedure for the analysis of the solid substances is as follows:

Aim: Qualitative analysis of given unknown Inorganic substance.

- (A) Physical properties
- (B) Dry tests for positive radicals or cations
- (C) Dry tests for negative radicals or anions
- (D) Conclusion from dry test:
 - (a) The probable positive radicals:
 - (b) The probable negative radicals :
- (E) Preparation of original solution
- (F) Wet tests for positive radicals or cations
- (G) Separation of cation into the groups (metals into groups)
- (H) Insoluble phosphate scheme (Only when substance is water insoluble)
- (I) Confirmative test for positive radicals or cations
- (J) Wet tests for negative radicals or anions (Only when substance is water soluble)

(K)	Confi	rmative test for neg	ative radicals	3					
(L)	Equat	tion for Cation and	Anion			ا ي∙			
(M)	Final	conclusion:	· ·	x	the c	250	FT . (o Jun.	•
	The g	iven unknown Inor	ganic substan	ce contains		1	e v efer	1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	
	(a)	Positive radical	:						
	(b)	Negative radical	:						

[A] PHYSICAL PROPERTIES:

NO.	TEST	OBSERVATION	INFERENCE
1	State	Crystalline (sandy powder)	All soluble substances like K ⁺ , NH ₄ ⁺ , Cl ⁻ , Br ⁻ , I ⁻ , NO ₃ ⁻ , SO ₄ ²⁻ etc. may be present.
		Amorphous (Soft powder)	:. Insoluble substances generally CO ₃ ²⁻ , S ²⁻ , PO ₄ ³⁻ , O ²⁻ etc. may be present.
2	Colour	White	Generally compounds of Al ³⁺ , Zn ²⁺ , Ca ²⁺ , Sr ²⁺ , Ba ²⁺ , Mg ²⁺ , K ⁺ , NH ₄ ⁺ etc. may be present.
		Black or Brown	MnO ₂ , CuO, Ni ₂ O, CdO, CuS, FeS, FeO, Fe ₂ O ₃ and some sulphides like SnS, FeS etc. may be present.
	5 pc	Green	:. Cr_{1}^{3+} , Fe^{2+} (ous), Ni^{2+} , $CuCl_2$, $FeSO_4$ may be present.
		Blue	Cu ²⁺ , Co ₃ (PO ₄) ₂ may be present.
		Yellow	:. CrO ₄ ⁻² , CdS, PbO, Bi ₂ O ₃ , sulphides of Cd ²⁺ , FePO ₄ , FeCl ₃ (Reddish yellow) may be present.
		Orange	∴ Cr ₂ O ₇ ²⁻ , Sb ₂ O ₃ may be present.
		Light pink	:. MnCl ₂ or MnSO ₄ may be present.
		Red (Dark Pink)	:. Co ²⁺ (hydrated), Pb ₃ O ₄ , Sb ₂ S ₃ etc. may be present.
		Reddish Brown	∴ Fe ₂ CO ₃ may be present.
3	Odour (smell)	Ammonical smell	∴ (NH ₄) ₂ CO ₃ , NH ₄ OH, Salts of NH ₄ ⁺ may be present.
	,	Rotten eggs (H ₂ S)	Soluble sulphides like K ₂ S, (NH ₄) ₂ S and insoluble sulphodes like ZnS, CuS etc. may be present.
4	Solubility in water	Soluble	All soluble substances like K ⁺ , NH ₄ ⁺ , Cl ⁻ , Br ⁻ , l ⁻ , NO ₃ , SO ₄ ²⁻ etc. may be present.
	Substance + H ₂ O (shake well) Heat if necessary	Insoluble	:. Insoluble substance like CO ₃ ²⁻ , S ²⁻ , PO ₄ ³⁻ , O ²⁻ etc. may be present.

				base or free acid may be			
5	Action on	Acidic	::	Salt of strong acid & weak base or free acid may be			
	Litmus paper	(Blue to red)		present.			
	nt.	Basic.	::	Salt of strong base & weak acid or free alkali may be			
		(Red to blue)		present.			
		Neutral	:	Salt of strong acid & strong base or			
		(no effect)		Salt of weak acid & weak base may be present.			
minn	Thi phy mich has a community of the comm						

		(no eff	ect)	Salt of weak acid & weak base ma		weak base may be present.
[B] D	RY TEST FOR POS	ITIVE RAD	ICAL:	1		
1	Heating in a dry tes [N.B.] (Take substance in	a dry		tance melts		Salt of K ⁺ , some salts containing water of crystallization may be present. Some nitrates like Pb(NO ₃) ₂ &
	test tube and heat it.) Test tube must be dibefore starting he	ry and eating,	-	duce cracking noise crepitates)		crystalline salts like KBr, KI may be present.
	bring wet litmus p starch paper, dichr paper, charred man etc.	omate part	Water condensing on the cooler part of the test tube		÷	Crystalline salt or some hydroxide may be present or hydroscopic like MgCl ₂ may be present.
20		cool	Sublimation takes place on the cooler part of the test tube. [Perform mirror test]			Hg ²⁺ , Sb ⁺³ , As ⁺³ , NH ₄ ⁺ etc. may be present.
	region of a release	Cha	nge of col	our	is an area	
2.		Blac	kening at l	high temperature		Cu ²⁺ , Mn ²⁺ , Ni ²⁺ , Co ²⁺ etc. may be present.
	, 13.00 			not and white when	:.	ZnO and many Zn ²⁺ salt may be present.
	./1	re la	Yellov	vish brown	::	Bi ₂ O ₂ and SnO ₂ may be present
	9. ≪ √		Y	ellow		PbO and some Pb ²⁺ salt may be present.
			В	rown	::-	Fe ³⁺ , CdO & Cd ²⁺ salts may be present.
				Green	∴	Cr^{3+} , CrO_4^{2-} , $Cr_2O_7^{2-}$ etc. may be present.
	white			e. white remains	.:	Al ³⁺ , Ca ²⁺ , Sr ²⁺ , Ba ²⁺ , Mg ²⁺ etc. may be present.
1	a a		ıtion of a			
		Evolu turns	ition of rec to wet star	ddish brown gas ch paper yellow	:.	Br may be present.

Evolution of reddish brown gas : NO ₃ may be present.	
no effect on wet starch paper	
Evolution of colourless or greenish yellow gas having pungent smell and bleaches (decolourized) wet litmus paper (Blue/red litmus paper turns to white)	a **
Evolution of NH ₃ gas recognized by smell, alkaline action on litmus, turmeric paper turns to brown and gives a white fumes with a drop of conc. HCl	
Evolution of Violet gas turns to wet starch paper blue :. I may be present.	
Evolution of a bad smelling gas (rotten eggs) turns to lead acetate paper black or silvery white S ²⁻ may be present: (N.B. Substance must be in water)	
White fumes of SO₂ gas acidic to litmus and having pungent smell ∴ SO₄²- (hydrated) may be	e present.
If the substance is yellow or orange and gives sparks on heating and becomes dark green on cooling.	present.
2 Mirror test White mirror ∴ Hg²+ may be present.	i.
(Perform when white sublimation take place) Orange sublimation ∴ Sb³+ may be present.	
Substance + Na ₂ CO ₃ + Garlic smell ∴ As ³⁺ may be present.	ed, P
Charcoal dust (Mix. thoroughly & heat the mixture in a dry test tube) Evolution of NH ₃ gas recognized by smell, alkaline action on litmus, turmeric paper turns to brown and gives a white fumes with a drop of conc. HCl	
3 Charcoal cavity test Make a cavity in the Deflagration ∴ NO₃⁻, Peroxide like I (Burning of charcoal) may be present.	BaO_2 , CrO_4^2
charcoal & moisten it with a drop of water. Put the mixture of Substance + Substance fuses and sinks in charcoal ∴ K ⁺ may be present.	, , , , , , , , , , , , , , , , , , ,
Na ₂ CO ₃ + Charcoal dust Sublimation ∴ NH ₄ ⁺ may be present.	
(mix well separately) in the cavity. Again put a drop of water on the mixture and heat it with White infusible mass which glows when heated etc. may be present. [PERFORM Co(NO ₃) ₂ TEST]	, Al ³⁺ , Mg ²⁴
blow pipe in the reducing Coloured infusible mass which :. Cu ²⁺ , Mn ²⁺ , Fe ²⁺ , Fe ³⁺	Co ²⁺ , Ni ²⁺ ,

	flame (yellow) for 8-10 minutes.	glows when heated [PERFORM BOREX BEAD TEST]	Cr^{3+} , CrO_4^{2-} , $Cr_2O_7^{2-}$ etc. may be present.
	r J. Tr.	Metallic Beads	(c) 00-1-1-
		White soft bead which marking on paper and yellow incrustation	∴ Pb ²⁺ may be present.
		Pinkish white bead with yellow incrustation	∴ Bi ²⁺ may be present.
	*	White bead not marking on paper	Sb ³⁺ may be present.
		Red scales	Cu ²⁺ may be present.
		Brown incrustation	:. Cd ²⁺ may be present.
		No bead, yellow when hot white when cold	∴ Zn ²⁺ may be present.
4	Cobalt nitrate test	Green mass	∴ Zn ²⁺ may be present.
	Put one or two drops of Co(NO ₃) ₂ solution on a		
	white infusible mass obtained in charcoal cavity and heat with blowpipe in	Blue mass	:. Al ³⁺ or PO ₄ ³⁻ may be present.
	reducing flame.	Pink mass	Mg ²⁺ may be present.
	[This test should be		:
	performe only when white infusible mass is obtained in charcoal cavity test]		
5	Borax Bead Test: Perform t	this test only when the substance is	coloured and leaves a coloured infusible

Borax Bead Test: Perform this test only when the substance is coloured and leaves a coloured infusible mass in charcoal cavity test.

(i.e. when any one of Cu^{2+} , Cr^{3+} , CrO_4^{2-} , $Cr_2O_7^{2-}$, Co^{2+} , Fe^{2+} , Fe^{3+} , Mn^{2+} , Ni^{2+} is suspected to be present).

Obser	vation	Inference		
Colour of the bead in the oxidizing flame	Colour of the bead in the reducing flame	Radicals		
dark yellow when hot Green when cold	Green (in hot & cold)	: Cr ³⁺ , CrO ₄ ²⁻ , Cr ₂ O ₇ ²⁻ etc. may be present.		
Blue (in hot & cold)	Blue (in hot & cold)	∴ Co ²⁺ may be present.		
Green when hot blue when cold	Colourless when hot Red or opaque when cold	∴ Cu ²⁺ may be present.		
Yellow	Light green	∴ Fe ²⁺ , Fe ³⁺ may be preset.		
Violet (Amethyst)	Colourless	∴ Mn ²⁺ may be present.		
Reddish brown when cold	Grey when cold	∴ Ni ²⁺ may be present.		

Flame Test: Omit this test if Cu²⁺, Pb²⁺, Sb³⁺, Bi³⁺ are detected from the preliminary dry test as the platinum wire is corroded by this radical. This test can however be performed by using a charred matchstick and when a green or bluish green flame is observed then the above metals are present.

	matchstick and when a green of bluish green hame is observed then the device metals are process.					
	Test	Obse	Observation			
	Prepare a paste of the substance in a minimum	Colour of the flame without cobalt glass.	Colour of the flame through the cobalt glass	Radicals		
_	quantity of conc. HCl	Brick red	Colourless	∴ Ca ²⁺ may be present.		
	and perform flame test in oxidizing flame (blue	Crimson	Crimson	∴ Sr ²⁺ may be present.		
	flame).	Light Green	Green	∴ Ba ²⁺ may be present.		
	,	Violet (Pink)	Crimson	∴ K ⁺ may be present.		
	i	Blue or greenish blue	-	∴ Cu ²⁺ may be present.		
		Bluish white (Pigeon Coloured)	7	∴ Pb ²⁺ , Sb ³⁺ , Bi ²⁺ etc. may be present.		
7	Test for NH ₄ ⁺ : (Sodium Hydroxide Test) Substance + NaOH (heat)	Evolution of NH ₃ ga alkaline action on litmo brown and gives a wh conc. HCl	NH ₁ ⁺ may be present.			

[C] DRY TEST FOR NEGATIVE RADICAL:

N.B. If given substance is insoluble in water then CO₃²⁻, S²⁻, O²⁻ and PO₄³⁻ (except salts of NH₄⁺ & K⁺) are to be present.

1	Test for CO ₃ ²⁻ and S ²⁻ : Substance + Dilute HCl (if substance is soluble in water and CO ₃ ²⁻ found to be present then K ⁺ or NH ₄ ⁺ is present)	Effervences of CO ₂ gas	: CO ₃ ²⁻ may be present.
	~		
2	Substance + Dilute HCl (heat)	Evolution of a bad smelling gas (rotten eggs) turns to lead acetate paper black or silvery white	∴ S ²⁻ may be present.
3	Test for Cl ⁻ , Br ⁻ , l ⁻ : Substance + MnO ₂ (powder) + Conc. H ₂ SO ₄ (heat) [N.B. Test tube must be dry]	Evolution of colourless or greenish yellow gas having pungent smell and bleaches (decolourized) wet litmus paper (Blue/red litmus paper turns to white)	∴ Cl ⁻ may be present.
	[N.D. Test tube must be dry]	Evolution of reddish brown gas turns to wet starch paper yellow	∴ Br may be present.
		Evolution of Violet gas turns to wet starch paper blue	∴ I may be present,

4	Test for NO ₃ :		
,3	(Perform when Br ⁻ or I ⁻ found to be absent) Substance + Conc. H ₂ SO ₄ + Cu foil (heat)	Evolution of Reddish brown gas & no effect on wet starch paper	:. NO ₃ may be present.
5	Test for PO ₄ ³⁻ : Substance + few drops of conc. HNO ₃ + (NH ₄) ₂ MoO ₄ in excess (Ammonium molybdate solution) warm slowly	Canary yellow precipitates	∴ PO ₄ ³⁻ may be present.
6	Test for CrO ₄ ²⁻ : If given Substance is yellow, add dil. HCl to its aqueous solution	Orange solution	∴ CrO ₄ ²⁻ may be present.
7	Test for Cr ₂ O ₇ ²⁻ : If given Substance is orange, add dilute NaOH solution to its aqueous solution	Yellow solution	∴ Cr ₂ O ₇ ²⁻ may be present.
8	Test for SO ₄ ⁻² : Crystalline substance + H ₂ O	Soluble	∴ SO ₄ ²⁻ may be present.
9	Test for O ·2:	All the above radicals are found to be absent and given substance is an amorphous and insoluble in water	∴ O ²⁻ may be present.

Note: Sulphates of Lead, Calcium, barium and strontium are sparingly soluble in water.

[D] CONCLUSION FROM DRY TEST:

[i]	Probable positive radical	:	
[ii]	Probable negative radical	:	

[E] PREPARATION OF ORIGINAL SOLUTION FOR WET TEST:

1	Substance + Distilled H ₂ O	Soluble	∴ O.S.
	(Shake well) (warm if it remains insoluble)		(Possibility of K ⁺ , NH ₄ ⁺ , NO ₃ ⁻ , SO ₄ ²⁻ , Cl ⁻ , Br , l ⁻) Or soluble PO ₄ ³⁺ solt like (NH ₄)HPO ₄ or K ₃ PO ₄)
2	Substance + dil. HCl	Soluble	$\therefore \text{ O.S.}$
	(heat if it is insoluble in cold)		(Possibility of CO ₃ ²⁻ , S ²⁻ , PO ₄ ³⁻)
3	Substance + Conc. HCl try to dissolve	Soluble	∴ O.S.
	as much of substance as possible, decant the solution and diluted with distill water.	v	(Possibility of CO ₃ ²⁻ , S ²⁻ , PO ₄ ³⁻)
	distill water.		·-

					1
4	Substance + HNO ₃ (heat) (prepared only when 1 st group is present)	Soluble	.:.	O.S. (Use for 1 st group radical)	

N.B.

- If insoluble sulphide (S^{2-}) is present (from dry test), boil the substance with conc. HCl for five minutes to expel H_2S gas. Then check with lead acetate paper till no effect on lead acetate paper i.e. all sulphide $(H_2S$ gas) is expel out.
- Lead acetate paper: Take a strip of filter paper and dip it into solution of lead acetate.
- If the powder (salt) is insoluble in water and PO_4^{3-} is detected then follow the "insoluble phosphate scheme" (Page no. 18) for the identification of the metallic radical after the 1^{st} and 2^{nd} groups.
- Use the solution (O.S.) prepared in acid for the determination of only positive radical.
- If O.S. is prepared in HCl, then 1st group is absent.
- If the substance is soluble in water and phosphate is detected then phosphate of K^+ , NH_4^+ may be present

[F] WET TEST FOR POSITIVE RADICALS:

Classification of positive radical into groups.

1	O.S.+ dilute HCl	White precipitates	 ∴ Ist group is present. i.e. Ag⁺, Hg⁺, Pb²⁺ is present.
2	O.S.+ dil. HCl + H ₂ S solution	Coloured precipitates	:. II nd group is present.
	(excess) + H ₂ O	Black or brown precipitates	Cu ²⁺ , Pb ²⁺ , Bi ³⁺ is present.
	(shake well)	Orange precipitates	∴ Sb ³⁺ is present.
		Yellow precipitates	∴ Cd ²⁺ is present.
	Test for PO ₄ ³⁻ : Perform phosp	phate test	Finite water y
		esent and substance is insoluble i e scheme (page no. 18) for rest of	in water and 1 st and 2 nd group are absent the elements.
	Test for PO ₄ ³⁻		, F :
	Substance + few drops of conc. HNO ₃ + (NH ₄) ₂ MoO ₄ in excess (Ammonium molybdate solution) warm	Canary yellow precipitates	∴ PO ₄ ³ is present.
	slowly.		
	[Now proceed as under in abs	ence of PO ₄ ³⁻]	
3	O.S. + NH ₄ Cl (excess in solid form) + 50% NH ₄ OH till solution became alkaline	White or coloured precipitates	∴ III rd A group is present. i.e. Al ³⁺ ,Fe ³⁺ ,Fe ²⁺ ,Cr ³⁺ present.
	and an analysis	White gelatinous sticky precipitates	∴ Al³+ is present.

		Reddish brown precipitates	∴ Fe ³⁺ (Ferric) is present.
		Green precipitates	∴ Fe ²⁺ (Ferrous) or Cr ³⁺ is present.
4	O.S. + NH ₄ Cl (excess in solid form) + 50% NH ₄ OH + H ₂ S solution (excess)	White, Buff or Black	i.e. Mn ²⁺ , Zn ²⁺ , Ni ²⁺ , Co ²⁺ present.
	solution (excess)	Buff precipitates Soluble in dilute acetic acid	∴ Mn ²⁺ is present.
		White precipitates Insoluble in dil. CH ₃ COOH	∴ Zn²+ is present.
		Black precipitates	∴ Ni ²⁺ or Co ²⁺ is present.
		O.S. is pink colour.	∴ Co ²⁺ is present.
		O.S. is green colour.	∴ Ni ²⁺ is present.
5	O.S. + NH ₄ Cl (Solid) + 50% NH ₄ OH + saturated solution of (NH ₄) ₂ CO ₃ excess.	White precipitates	∴ IV th Group is present. i.e. Ca ²⁺ ,Sr ²⁺ , Ba ²⁺ is present.
6	O.S. + NH ₄ Cl (solid) + 50% NH ₄ OH + Na ₂ HPO ₄	White precipitates	 ∴ Vth A Group is present. i.e. Mg²⁺ is present.
7	Test for V th B Group:	All the above groups (i.e. from group I to V th A) are found to be absent and also given substance is soluble in water	∴ V th B Group is present. i.e. K ⁺ , NH ₄ ⁺ is present.

[G] SEPARATION OF CATION INTO GROUPS:

GROUP I:

Separation of Ag+, Hg+ & Pb2+

1	O.S. + dil. HCl	White precipitates	∴ Ag ⁺ , Hg ⁺² & Pb ²⁺ is present.
2	Take small quantity of	XX 71	∴ Pb ²⁺ is present.
	white precipitates of above test + NH ₄ OH (excess)	3371.14	
		White precipitates Calubi	∴ Ag ⁺ is present.

C.T. for Pb⁺²

l O.S. + K ₂ CrO ₄	Yellow precipitates of (PbCrO ₄) soluble in hot NaOH	∴ Pb ⁺² is confirmed.
2 O.S. + KI solution	Yellow precipitates of (PbI ₂) soluble in hot (boiling) water, recrystallises into golden yellow needle (plates) on cooling	∴ Pb ⁺² is confirmed.

			-		
3	O.S. + NaOH solution	White precipitates of Pb(OH) ₂ soluble in excess of NaOH	÷	Pb ⁺² is confirmed.	
GRO	UP II :				
Separ	ration of Cu ²⁺ , Pb ²⁺ , Bi ³⁺ (Black	precipitation)		4 + 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	
1	O.S. + NaOH solution	Bluish precipitates turns to black on heating	:.	Cu ²⁺ is present.	
		White precipitates insoluble in excess of NaOH	:	Bi ³⁺ is present.	
		White precipitates soluble in excess of NaOH		Pb ²⁺ is present.	
2	O.S. + dilute H ₂ SO ₄	White precipitates,	••	Pb ²⁺ is present	
Separ	ration of As ³⁺ , Cd ²⁺ , Sb ³⁺ (Yello	w precipitation)	÷		
1	O.S. + dil. HCl + H ₂ S solution (excess) + H ₂ O	Yellow precipitates soluble in # conc. HCl (shake well)	••	Cd ²⁺ is present. (2.20 + 2.30)	
	(excess) shake well	Orange precipitation	:	Sb ⁺³ is present.	
C.T.	for Cu ²⁺			. Fig. 1 2 1 1	
1	O.S. + NaOH	Blue precipitates turns to black on heating insoluble in excess of NaOH	2.	Cu ²⁺ is confirmed.	
2	O.S. + NH₄OH	Bluish precipitates Soluble in excess of NH ₄ OH giving deep blue solution.	••	Cu ²⁺ is confirmed.	
3	O.S. + KI	Dirty green precipitates of (Cu ₂ I ₂) appearing brown due to liberation of I ₂	••	Cu ²⁺ is confirmed.	
4	O.S. + Potassium Ferrocyanide [K ₄ Fe(CN) ₆]	Reddish brown precipitates soluble in NH ₄ OH and insoluble in dilute acid	l .	Cu ²⁺ is confirmed.	
C.T.	for Bi ³⁺				
1	O.S. + NaOH	White precipitates of [Bi(OH) ₃] insoluble in excess of NaOH		∴ Bi ³⁺ is confirmed.	
2	O.S. + KI	Brown precipitates of (Bil ₃) turns to orange solution on addition of excess of KI			
3	O.S. + H ₂ S solution excess	Brown precipitates of (Bi ₂ S ₃) insoluble in tartaric acid (cold di	lute)	∴ Bi ³⁺ is confirmed.	

C.T. for Pb⁺²

1	O.S. + K ₂ CrO ₄	Yellow precipitates of (PbCrO ₄) soluble in hot NaOH	∴ Pb ⁺² is confirmed.
2	O.S. + - KI	Yellow precipitates of (PbI ₂) soluble in hot water, recrystallises into golden yellow needle on cooling	∴ Pb ⁺² is confirmed.
3	O.S. + NaOH	White precipitates of Pb(OH) ₂ soluble in excess of NaOH	∴ Pb ⁺² is confirmed.
4	O.S. + dil. H ₂ SO ₄	White precipitates	∴ Pb ⁺² is confirmed.

C.T. for Cd²⁺

1	O.S. + Na ₂ CO ₃	White precipitates of CdCO ₃	Cd ²⁺ is confirmed.
		s · _ ; * ;	rita i i i i i i i i i i i i i i i i i i
2	O.S. + NaOH	White precipitates of	∴ Cd ²⁺ is confirmed.
,		[Cd(OH) ₂] insoluble in excess	*
	ys with my the	of NaOH	and the
3	O.S. + HgCI ₂	White precipitates of (Hg ₂ CI ₂)	∴ Cd ²⁺ is confirmed.
	10 0 0 0 0	turns to gray in excess of	
		HgCI ₂	1 10

C.T. for Sb³⁺

1	O.S. + H ₂ O	White precipitates of (SbOCl) soluble in tartaric acid	∴ Sb ³⁺ is confirmed.
2	O.S. + AgNO ₃ + NaOH (excess)	Black precipitates insoluble in ammonia	∴ Sb³+ is confirmed.
3	O.S. + NaOH	White precipitates of Sb ₂ O ₃ soluble in excess of NaOH	∴ Sb ³⁺ is confirmed.

GROUP III A:

Separation of Al^{3+} , Fe^{3+} , Fe^{2+} and Cr^{3+}

1	O.S. + NH_4Cl + 50% NH_4OH (excess)		White gelatinous precipitates (O.S. is colourless)	∴ Al ³⁺ is present.
	, xv 1	1 1 1	Reddish brown precipitates (O.S. is yellow)	∴ Fe ³⁺ (Ferric) is present.
	• + ,9 - , - +)	i a	Dirty green precipitates turns to brown on surface due to exposure in air	∴ Fe ²⁺ (Ferrous) is present.
	. 2 .	73	Bluish green precipitates (O.S. is bluish green)	∴ Cr ³⁺ is present.

Separation of Fe^{2+} from Cr^{3+}

O.S. + 2-3 drops of Conc.	Reddish brown precipitates	∴ Fe ²⁺ is present.
HNO ₃ boil and cool + 50 %	Bluish green precipitates	∴ Cr ³⁺ is present.
NH₄OH (excess)	Yel cipit (PE	

C.T. for Al³⁺

1	O.S. + NaOH (shake well)	White gelatinous sticky precipitates soluble in excess of NaOH and reprecipitated by NH ₄ Cl	∴ Al ³⁺ is confirmed.
2	O.S. + NH ₄ OH + H ₂ S solution	White precipitates of Al(OH) ₃ soluble in excess	∴ Al ³⁺ is confirmed.
3	O.S. + K ₄ Fe(CN) ₆	No precipitates	∴ Al ³⁺ is confirmed.

C.T. for Fe²⁺

1	O.S. + K ₄ [Fe(CN) ₆] (Potassium ferrocyanide)	Bluish white precipitates of Fe ₂ [Fe(CN) ₆]	∴ Fe ²⁺ is confirmed.
2	O.S. + K ₃ [Fe(CN) ₆] (Potassium ferricyanide)	Dark blue precipitates of Fe ₃ [Fe(CN) ₆] ₂	∴ Fe ²⁺ is confirmed.
3	O.S. + NaOH or NH ₄ OH	Dirty green precipitates of Fe(OH) ₂ turns to reddish brown on exposure in air	∴ Fe ²⁺ is confirmed.
4	O.S. + Dilute H ₂ SO ₄ + dilute KMnO ₄ (Two drops) (shake well)	Colour of KMnO ₄ discharged (disappear)	∴ Fe ²⁺ is confirmed.

C.T. for Fe³⁺

1	O.S. + K ₄ Fe(CN) ₆ (Potassium ferrocyanide)	Deep blue precipitates of Fe ₄ [Fe(CN) ₆] ₃	*	Fe ³⁺ is confirmed.
2	O.S. + K ₃ Fe(CN) ₆ (Potassium ferricyanide)	Brown colour solution	•	Fe ³⁺ is confirmed.
3	O.S. + KCNS	Blood red colour K ₃ [Fe(SCN) ₆]	••	Fe ³⁺ is confirmed.
4	O.S. + Na ₂ HPO ₄	Yellowish white precipitates of FePO ₄ insoluble in dilute acetic acid.	:.	Fe ³⁺ is confirmed.
5	O.S. + NaOH or NH ₄ OH	Reddish brown precipitates of Fe(OH) ₃ insoluble in excess	:	Fe ³⁺ is confirmed.

C.T. for Cr³⁺

1	O.S. + PbO ₂ + NaOH (boil) + acetic acid	Yellow precipitates of PbCrO ₄	∴ Cr ³⁺ is confirmed.
2	O.S. + NH₄OH	Bluish green precipitates of Cr(OH) ₃ soluble in excess of NH ₄ OH	∴ Cr ³⁺ is confirmed.
3	O.S. + NaOH	Bluish green precipitates of Cr(OH) ₃	∴ Cr ³⁺ is confirmed.

GROUP III B:

Separation of Zn2+, Mn2+, Ni2+, Co2+

l	O.S. + NH ₄ Cl (excess)	White precipitates insoluble in acetic acid	:.	Zn ²⁺ is present.
	+ 50% NH₄OH + H₂S solution (excess)	Pinkish or white or buff coloured precipitates soluble in acetic acid	÷	Mn ²⁺ is present.
		Black precipitates	÷	Ni ²⁺ or Co ²⁺ is present.
2	O.S. + NaOH	White precipitates soluble in excess of NaOH	••	Zn ²⁺ is present.
	(shake well)	White precipitates insoluble in excess NaOH	÷	Mn ²⁺ is present.
	19. 10 2	Light green precipitates insoluble in excess of NaOH (O.S. must be green)	÷	Ni ²⁺ is present.
		Blue precipitates turns to pink on heating	::	Co ²⁺ is present.
		(O.S. must be pink)	,	ay. The same

C.T. for Zn²⁺

1	O.S. + K ₄ Fe(CN) ₆ (Potassium ferrocyanide)	White precipitates of Zn ₂ [Fe(CN) ₆] insoluble in dil. HCl	∴ Zn ²⁺ is confirmed.
2	O.S. + NaOH or NH ₄ OH	White precipitates of Zn(OH) ₂ soluble in excess	∴ Zn ²⁺ is confirmed.
3	O.S. + K ₃ Fe(CN) ₆ (Potassium ferricyanide)	Orange precipitates soluble in dil. HCl	∴ Zn ²⁺ is confirmed.

C.T. for Mn²⁺

1	O.S. + K ₄ Fe(CN) ₆ (Potassium ferrocyanide)	Pinkish white precipitates of Mn ₂ [Fe(CN) ₆] soluble in dil. HCl (no precipitation if O.S. is acidic)	Mn ²⁺ is confirmed.
2	O.S. + PbO ₂ + conc. HNO ₃ boil, cool and dilute with water and allow it to settle.	Pink colour of KMnO ₄	∴ Mn ²⁺ is confirmed.

C.T. for Co²⁺

	O.S. + sodium nitrite (NaNO ₂) + CH ₃ COOH + KI	Yellow precipitates of [K ₃ (Co(NO ₂) ₆)]	∴ Co ²⁺ is confirmed.
2	O.S. + Na ₂ HPO ₄	Violet precipitates of [Co ₃ (Po ₄) ₂]	∴ Co ²⁺ is confirmed.
3	O.S. + NaOH	Blue precipitates of [Co(OH) ₂] insoluble in excess of NaOH & turns to pink on heating	∴ Co ²⁺ is confirmed.

C.T. for Ni²⁺

1	O.S. + Br ₂ water + NaOH in excess	Black precipitates of (Ni ₂ O ₄ .3H ₂ O)	∴ Ni ²⁺ is confirmed.
2	O.S. + Dimethyl glyoxime (C ₄ H ₇ C ₂ N ₂) + 50% NH ₄ OH excess (shake well)	Deep red precipitates	:. Ni ²⁺ is confirmed.
3	O.S. + NH ₄ OH (2-3 drops)	Light green precipitates soluble in excess of NH ₄ OH	Ni ²⁺ is confirmed.

GROUP IV:

Separation of Ca²⁺, Ba²⁺ & Sr²⁺

1	O.S. + CaSO ₄ solution	White precipitates immediately	∴	Ba ²⁺ is present.
	(O.S. must be concentrate)	White precipitates formed on warming or appearing slowly on standing	÷	Sr ²⁺ is present.
		No precipitates	··	Ca ²⁺ is present.

C.T. for Ca²⁺

1	$O.S. + K_2CrO_4$	No precipitates	∴ Ca ²⁺ is confirmed.
2	O.S. + Ammonium oxalate [(NH ₄) ₂ C ₂ O ₄]	White precipitates of CaC ₂ O ₄ insoluble in CH ₃ COOH	∴ Ca ²⁺ is confirmed.
3	O.S. + dilute H ₂ SO ₄	White precipitates soluble in ammonium sulphate	∴ Ca ²⁺ is confirmed.

C.T. for Ba²⁺

1	O.S. + K ₂ CrO ₄	Yellow precipitates of (BaCrO ₄) insoluble in dilute CH ₃ COOH but soluble in dil. HCl	∴ Ba ²⁺ is confirmed.
2	O.S. + Ammonium oxalate [(NH ₄) ₂ C ₂ O ₄]	White precipitates of BaC ₂ O ₄ soluble in CH ₃ COOH on heating	∴ Ba ²⁺ is confirmed.
3	O.S. + dil. H ₂ SO ₄	White precipitates soluble in all acids like HCl, HNO ₃ etc.	∴ Ba ²⁺ is confirmed.

C.T. for Sr²⁺

1	O.S. + K ₂ CrO ₄	Yellow precipitates of (SrCrO ₄) on heating	∴ Sr ²⁺ is confirmed.
	41.5	soluble in CH ₃ COOH	
2	O.S. + $(NH_4)_2C_2O_4$	White precipitates of (SrC ₂ O ₄) soluble in CH ₃ COOH on heating	Sr ²⁺ is confirmed.
3	O.S. + dil. H ₂ SO ₄	White precipitates of (SrC ₂ O ₄) after few minutes	∴ Sr ²⁺ is confirmed.

GROUP VA:

Identification of Mg⁺²

1	O.S. + NH ₄ Cl + NH ₄ OH +	White precipitates	::	Mg ²⁺ is present.	
	Na ₂ HPO ₄				

C.T. for Mg²⁺

1	HYPOIODIDE TEST:		
	O.S. + KI + NaOH + I ₂ solution till solution becomes yellow	Reddish brown precipitates of [Mg(IO) ₂]	∴ Mg ²⁺ is confirmed.
2	O.S. + (NH ₄) ₂ C ₂ O ₄ (Ammonium oxalate)	White precipitates of (MgC ₂ O ₄) soluble in CH ₃ COOH	∴ Mg ²⁺ is confirmed.
3	O.S. + NH ₄ OH	White precipitates soluble in NH ₄ Cl and HCl	Mg ²⁺ is confirmed.

Hypoiodide solution: To NaOH solution add equal amount of KI and few drops of l_2 solution till yellow colour is obtained.

If no precipitation obtained in all above groups (i.e. from group I^{st} to $V^{th}A$) then V^{th} B group (K^+, NH_4^+) is consider to be present.

GROUP Vth B:

Separation of K+ and NH4+

1	O.S. + Freshly prepared Sodium cobaltinitrite	Yellow precipitates	:.	K ⁺ or NH ₄ ⁺ is present.
2	O.S. + NaOH (heat)	Evolution of NH ₃ gas recognized by smell, alkaline action on litmus, turmeric paper turns to brown and gives a white fumes with a drop of conc. HCl		NH ₄ ⁺ is present.
		No evolution of NH ₃ gas	·:	K ⁺ is present.
3	O.S. + Nessler's reagent	Reduisii blowli precipitates		NH ₄ ⁺ is present.

Sodium cobaltinitrite { $Na_3[Co(NO_2)_6]$ }: Take small quantity of solid sodium nitrite ($NaNO_2$), add few drops of cobaltnitrate[$Co(NO_3)_2$] solution and shake it. Add dilute acetic acid solution.

Now check this sodium cobaltinitrite solution with NH₄Cl solution: Take one drop of NH₄Cl solution add sodium cobaltinitrite solution in excess. It gives yellow precipitates, means sodium cobiltinitrite solution is true. If no precipitation means that prepared sodium cobaltinitrate solution is not true.

Nessler's reagent { $K_2[HgI_4]$ }: Take one drops of $HgCl_2$ solution, add to it KI solution drop wise with shaking till precipitates just dissolve. Add an equal volume of NaOH and shake it.

Now check with NH_4Cl solution. (one drop of NH_4Cl + nessler's reagent (excess) give yellowish brown precipitates, If no precipitation means, prepared nessler's reagent is not true.)

N.B. If the substance is soluble in water and phosphate or carbonate is found to be present. Then phosphate or carbonate of K^+ and NH_4^+ may be present. All other phosphate and carbonate are insoluble in water.

C. T. for K+

1	O.S. + Freshly prepared Sodium cobaltinitrite	Yellow precipitates	∴ K ⁺ is confirmed.
2	O.S. + Picric acid	Yellow precipitates	∴ K ⁺ is confirmed.
	(shake well)	81	
3	O.S. + Tartaric acid	White precipitates	∴ K ⁺ is confirmed.
	(shake well)		
4	O.S. + Nessler's reagent	No precipitates	∴ K ⁺ is confirmed.

C. T. for NH₄⁺

1	O.S. + Nessler's reagent	Reddish-brown precipitates	\cdot	NH ₄ ⁺ is confirmed.
2	O.S. + NaOH (heat)	Evolution of NH ₃ gas recognized by smell, alkaline action on litmus, turmeric paper turns to brown and gives a white fumes with a drop of conc. HCl		NH ₄ ⁺ is confirmed.

[H] INSOLUBLE PHOSPHATE SCHEME:

If given substance is insoluble in water and 1st and 2nd group are found to be absent and PO₄-3 is present then follow the INSOLUBLE PHOSPHATE SCHEME as given below.

Identification of positive radicals of group III^{rd} A, III^{rd} B, IV^{th} and V^{th} group

1	O.S. + NH ₄ Cl (solid) + 50%	Black precipitates	::	Fe ²⁺ , Fe ³⁺ , Ni ²⁺ & Co ²⁺ is
	NH ₄ OH + H ₂ S solution in	ac I		present.
	excess	White precipitates		Al ³⁺ , Zn ²⁺ , Ca ²⁺ , Ba ²⁺ , Sr ²⁺ & Mg ²⁺ is present.
		Green precipitates	::	Cr ³⁺ is present.
		(O.S. is bluish green colour)		
	ÿ.	Buff coloured precipitates soluble in acetic acid		Mn ²⁺ is present.

Separation of black precipitation of Fe2+, Fe3+, Ni2+, Co2+

1	O.S. + NaOH sol ⁿ . till just alkaline (excess)	Dirty green precipitates changing to reddish brown on exposure of air and insoluble in excess of NaOH	∴ Fe ²⁺ is present.
		Reddish brown precipitates	∴ Fe³+ is present.
	de dan wan god hein	Light green precipitates (O.S. is green)	∴ Ni ²⁺ is present.
	mercy after reflorant to	Bluish precipitates (O.S. is red)	∴ Co²+ is present.

Separation of white precipitation of Al3+, Zn2+, Ca2+, Ba2+, Sr2+ & Mg2+

1	O.S. + 2 – 3 drops of NaOH solution	(1) White precipitates soluble in excess of NaOH	:. Al ³⁺ or Zn ²⁺ is present.
	K* is combinated	(a) Reprecipitated on addition of NH ₄ OH solution	:. Al ³⁺ is present.
-	Lexquinozzi A (1)	(b) No reprecipitation on addition of NH ₄ OH solution	∴ Zn ²⁺ is present.
	E is confirmed a Confirmed a Confirmed	(2) White precipitates insoluble in excess of NaOH	.: Ca ²⁺ , Ba ²⁺ Mg ²⁺ & Sr ²⁺ is present.

Identification of Ca2+, Ba2+ Mg2+ & Sr2+

1	O.S. + CaSO ₄ solution	White precipitates immediately	∴ Ba ²⁺ is present.
	no paper of a Mark to continue	White precipitates slowly or on warming	∴ Sr ²⁺ is present.
	planting	No precipitation even on warming	∴ Ca ²⁺ , Mg ²⁺ is present.
2	O.S. + Na ₂ CO ₃ till	Reprecipitates	∴ Ca ²⁺ is present.
	precipitates formed + Acetic acid till just dissolves + ammonia oxalate (excess)	No reprecipitation	∴ Mg ²⁺ is present.

[I] CONFIRMATIVE TEST FOR POSITIVE RADICAL:

Write confirmative test from earlier pages, whatever is applicable.

[J] WET TEST FOR NEGATIVE RADICALS (ANION):

- N.B: (a) When the substance is soluble in water, prepare its original solution in distilled water and detect the negative radical as under.
 - (b) When given substance is insoluble in water then negative radical can be detected only from dry test.

1	O.S. + AgNO ₃	(A) White or yellowish white precipitates	∴ Cl ⁻ , Br ⁻ , l ⁻ is present.
1	(shake well)	insoluble in dil. HNO ₃	
	(snake well)	Take small quantity of above white	
		precipitates of (A) + NH ₄ OH	i akas ,cens
	20	(shake well)	
	,	(i) White precipitates soluble in NH₄OH	∴ Cl ⁻ is present.
		(ii) Yellowish precipitates soluble in excess of NH ₄ OH	∴ Br is present.
		(iii) Yellowish precipitates insoluble in excess of NH ₄ OH	∴ I is present.
		(B) Yellowish precipitates soluble in dilute HNO ₃ without any effervences	∴ PO ₄ ³⁻ is present.
		(C) Whitish grey precipitates soluble in dilute HNO ₃ with effervences of CO ₂ gas	∴ CO ₃ ²⁻ is present.
		(D) Scarlet red precipitates soluble in dil. HNO ₃	:. CrO ₄ ²⁻ or Cr ₂ O ₇ ²⁻ is present.
		(i) O.S. is yellow, turns orange on addition of dilute HCl	∴ CrO ₄ ²⁻ is present.
		(ii) O.S. is orange, turns yellow on adding NaOH	\therefore $\operatorname{Cr}_2\operatorname{O}_7^{2-}$ is present.
		(E) No precipitates	∴ SO ₄ ²⁻ or NO ₃ is present.
2	O.S. + BaCl ₂	White precipitates insoluble in dilute HCl	∴ SO ₄ ²⁻ is present.
	solution	No precipitates	∴ NO ₃ is present.

N. B.: (a) Use $Ba(NO_3)_2$ instead of $BaCl_2$ if $\ I^{st}$ group is present.

(b) Wet tests and C.T. for (-Ve) negative radical cannot be performed when the given substance is insoluble in water and it can be concluded only from the dry test.
[Insoluble CO₃²⁻, O²⁻, PO₄³⁻ and S⁻² can be detected from dry test.]

1	Substance + dilute HCl	Effervences of CO ₂ gas	∴ CO ₃ ²⁻ is confirmed.
2	Substance + Conc. HCl (heat)	Evolution of a bad smelling gas (rotten eggs) turns to lead acetate paper black or silvery white	∴ S ²⁻ is confirmed.
3	Substance + Conc. HNO ₃ + Ammonium Molybdate (warm it)	Canary yellow precipitates	∴ PQ ₄ ³⁻ is confirmed.

4	Test for O ⁻² :	All the above negative radicals are	∴ O ²⁻ is confirmed.
		found to be absent and the given	(shake well) inso
		substance is an amorphous and	let.
	-	insoluble in water	
	CO ₃ the confirmed.	109°C	O S Lead aceta
	CONFIRMATORY TEST FO	R NEGATIVE RADICALS:	an and an and an and an
1	O.S. + Lead acetate	White precipitates soluble on	∴ Cl is confirmed.
		heating and reprecipitated on cooling	and August
		300	1377
2	Substance+ MnO ₂ (Powder)	Evolution of greenish yellow gas	∴ Cl ⁻ is confirmed.
	+ Conc. H ₂ SO ₄ and heat	having pungent smell, bleaches the	HM.
		litmus paper (Blue/red litmus paper	
	manufactor a re-	turns to white)	(a) Yes
C. T.	for Br	audit of stations authorized uses detail	AME ON
1	O.S. + Lead acetate	White precipitates soluble on	∴ Br is confirmed.
		heating and reprecipitated. on cooling	
2	O.S. + MnO ₂ + Con. H ₂ SO ₄	1/ 1 The se state to see a school beautiful	
2	in present	Reddish brown gas turns to wet	∴ Br is confirmed.
	(Heat)	starch paper yellow (orange)	116
3	O.S. + CHCl ₃ + Cl ₂ water	Yellow globules at bottom of test	Br is confirmed.
	(Shake well)	tube	Bi is committed.
С. Т.	for I	ashiriqi oqiq	on (3)
	00.1.1	Ly will supul. 1. J follow . A supul since	D 000000 200 200
1	O.S. + Lead acetate	Yellow precipitates soluble on	∴ Γ is confirmed.
	NO is present.	heating and becomes golden yellow	normicz
2	$O.S. + MnO_2 + Con. H_2SO_4$	Violet gas turns to wet starch paper	∴ I is confirmed.
	(heat)	blue	I is commined.
2	4 contribution given substance 4	Wielet elebula and h	01 07 09 100 mm 2 mm 2 mm
3	O.S. + $CHCl_3 + Cl_2$ water	Violet globules at the bottom of the	∴ I is confirmed.
	(Shake well)	test tube	30 300 4
4	O.S. + HgCl ₂	Deep red precipitates soluble in	· Lison C
	* Dearthmeet (O)	excess of HgCl ₂	: I is confirmed.
C T	6 po -3		
C. I.	for PO ₄ ⁻³		
1	O.S. + Conc. HNO ₃ +	Canary yellow precipitates	∴ PO ₄ ³⁻ is confirmed.
	NH ₄) ₂ MoO ₄ (Ammonium	MAN year (Control of the Control of	is confirmed.
	molybdate) in excess (warm i	t)	
2	O.S. + Lead acetate	White precipitates	OHI-
	o.s. 1 Lead acctate	winte precipitates	∴ PO ₄ ³⁻ is confirmed.
3	O.S. + Magnesia mixture.	White precipitates	
	THE RESIDENCE OF THE PARTY OF T	THE RESERVE TO BE A PERSON OF THE PERSON OF	∴ PO ₄ ³⁻ is confirmed.

C. T. for CO_3^{-2}

1	O.S. + BaCl ₂ or CaCl ₂	White precipitates soluble in acids with effervences of CO ₂ gas	∴ CO ₃ ²⁻ is confirmed.
 2	O.S. + Lead acetate	White precipitates	∴ CO ₃ ²⁻ is confirmed.
3	O.S. + dilute HCl	Effervences of CO ₂ gas turns to	\therefore CO ₃ ² is confirmed.
	(Shake well)	lime water milky	·

C. T. for NO₃

1	O.S. + Cu foil + Conc. H ₂ SO ₄ heat slowly	Reddish brown fumes of NO ₂ gas and no effect on wet starch paper	∴ NO ₃ is confirmed.
2	O.S. + Zn dust + NaOH (boil)	Evolution of NH ₃ gas turns to wet red litmus paper blue	∴ NO ₃ is confirmed.
3	Ring test: O.S. + Conc. H ₂ SO ₄ cool and add Freshly prepared FeSO ₄ solution from the side of the test tube.	Brown ring at the junction of two liquids	∴ NO ₃ is confirmed.

C. T. for SO₄⁻²

1	O.S. + BaCl ₂		\therefore SO ₄ ²⁻ is confirmed.
	,	acid	
2	O.S. + Lead acetate	White precipitates soluble in ammonium acetate solution	∴ SO ₄ ²⁻ is confirmed.
3	O.S. + Pb(NO ₃) ₂	White precipitates insoluble in acid and soluble in hot ammonium acetate (NH ₄ Ac)	∴ SO ₄ ²⁻ is confirmed.

C. T. for CrO₄⁻²

1	Colour of original solution	Yellow	∴ CrO ₄ ²⁻ is confirmed.
2	O.S. + BaCl ₂ solution	Yellow precipitates soluble in dilute HNO ₃	∴ CrO ₄ ²⁻ is confirmed.
3	O.S. + Dilute HCl	Orange colour	CrO ₄ ²⁻ is confirmed.
4	O.S. + Lead acetate [Pb(Ac) ₂]	Yellow precipitates	∴ CrO ₄ ²⁻ is confirmed.

C. T. for $Cr_2O_7^{-2}$

1	Colo	our of Orig	inal solution	Orange				<i>:</i> .	\therefore Cr ₂ O ₇ ²⁻ is confirmed.		
2	O.S	. + BaCl ₂		Yellow precipitates soluble in dilute :. $Cr_2O_7^{2-}$ is confirmed.							
		*17)ffCr *	HNO ₃	، بند	211 -	Ü. I			.K ,	
3	O.S.	+ Lead ac	etate [Pb(Ac) ₂]	Yellow precipitates			<i>:</i> .	$Cr_2O_7^{2-}$ is c	onfirmed.		

[L] EQUATION FOR CATION:

[L] EQUATION FOR ANION:

[M] FINAL CONCLUSION:

The given unknown substance contains:

Positive radical

Negative radical

Formula of the substance:

Name of the substance

PREPARATION OF THE REAGENTS

Sodium cobaltinitrate reagent Na ₃ [Co(No ₂) ₆]	Take equal volume of Co(NO ₃) and NaNO ₂ solution. Add dilute acetic acid to it. [Test with KI solution. It will give yellow precipitates]
Sodium Hypo iodide reagent [NaOI]	To NaOH solution add an equal volume of Kl solution and then add few drops of l_2 solution till permanent yellow colour is obtained.
Nessler's reagent K ₂ [HgI ₄]	Take few drops of HgCl ₂ and add KI solution to it till precipitates just dissolved. And then add an equal volume of NaOH solution.
Magnesia mixture	Take 5 gm of MgCl ₂ and dissolve in water (4-5 ml) to this solution add NH ₄ OH till it become alkaline. Then add NH ₄ Cl solution till the precipitates just dissolve.
Cobalt acetate [Co(CH ₃ Coo) ₂]	Take cobaltinitrate solution [Co(NO ₃) ₂] and add acetic acid to it.
Neutral FeCl ₃	Take FeCl ₃ solution and add NH ₄ OH to it, till reddish brown precipitates just form to start and then add water.

SOME CATION AND ANION CAN BE DIRECTLY CONFIRMED BY USING FOLLOWING TEST

- Ni²⁺: Substance + 50 % NH₄OH excess + DMG shake well, it gives scarlet red precipitates (Generally given substance is green or bluish green)
- Cu²⁺: Substance + 50 % NH₄OH shake well, it gives deep blue colouration. (Generally given substance is blue)
- Fe²⁺: Substance + Conc. NHO₃ + KCNS solution shake well, it gives blood red colour.
- Mg²⁺: Substance + Hypoiodide solution shake well, gives reddish brown precipitation.
- PO₄⁻³: Substance + Ammonium molybledum + conc.HNO₃ give cannery yellow precipitates.

 [If substance is blue (copper phosphate) then yellow precipitation is not clear yellow.]
- I -: Substance + conc. H₂SO₄ heat gently gives bluish violet gas.
- NO₃⁻: Substance + KI (solid) + starch solution + Zn dust gives blue colour.
- S⁻²: Substance (Insoluble in water) + Conc. HCl, heat, evolution of a bad smelling gas (rotten eggs) turns to lead acetate paper black or silvery white.
- SO_4^{-2} : Substance + BaCl₂ gives white precipitates...